

Self-Seal Sterilization Pouches

Sterilisation pouches with a self-adhesive strip to allow closure without the need for heat sealing equipment. These pouches are available in a variety of sizes and are printed with batch information as well as indicators for steam and EO processing. They come with a CE certificate and are manufactured in a class 10,000 clean room controlled to ISO 1384:2003.



Product Features

- Heavyweight 70gsm paper
- Blue tinted transparent 2-ply laminate film allows easy viewing of the contents of the pouch
- Suitable for both steam and ethylene oxide sterilisation
- Each pouch printed with Class A process indicators for both steam and ethylene oxide sterilisation
- Each pouch printed with a lot number for full manufacturing traceability
- Self Seal strip eliminates the need for heat seal equipment and provides a validated seal
- Off-Set webs for easy opening
- Tack seals to prevent curling and reduce external contamination
- Comprehensive size range

Available Size Options:-

Code	Nominal Size	Box Qty	Case Quantity
11110	2¼" x 4" / 60mm x 100mm	Box 200	10 boxes of 200 per case
11120	3½" x 5½" / 90mm x 140mm	Box 200	10 boxes of 200 per case
11128	3½" x 8" / 90mm x 200mm	Box 200	10 boxes of 200 per case
11130	3½" x 9" / 90mm x 230mm	Box 200	10 boxes of 200 per case
11140	5¼" x 10" / 135mm x 255mm	Box 200	10 boxes of 200 per case
11144	5¼" x 14" / 135mm x 360mm	Box 200	10 boxes of 200 per case
11150	7½" x 12" / 190mm x 305mm	Box 200	5 boxes of 200 per case
11153	7½" x 13" / 190mm x 330mm	Box 200	5 boxes of 200 per case
11185	12" x 15" / 305mm x 381mm	Box 200	5 boxes of 200 per case

We are also able to manufacture custom sizes and undertake private label production, by arrangement.

Minimum order quantity – 3 boxes.