


## X Ray Analysis of Medical Devices

Non destructive, non contact X-ray analysis of medical devices is now available. Medical Engineering Technologies had signed an exclusive agreement with Scan Vision for UK and Ireland distribution of the Rayonlab: High resolution X-Ray Station for Plastic Components and Packages.

This CE marked inherently safe unit delivers X-ray photographs showing hidden details of products and inside packaging. A robot (4 axis) positions the device for examination within the test chamber and can move it or rotate it for video footage. The radiation used is low energy and is contained within the chamber with full safety interlocks. Single phase power is required.

This system is ideal for analysis of: guide-wires and stents, product hidden within packaging, the internal components of drug delivery devices, the mating of surfaces and for detecting metallic contamination.



Inhaler assembly

Medical Engineering Technologies Ltd  
Yew Tree Studios, Stone Street, Stanford North, Ashford, Kent TN25 6DH, UK.

**Tel.:** +44 (0)8454 588924

**Fax:** +44 (0)1303 813592

**E-mail:** [solutions@met.uk.com](mailto:solutions@met.uk.com)

**Web:** [www.met.uk.com](http://www.met.uk.com)

For immediate release: Ashford, 7<sup>th</sup> April 2011